BI SoberTech[™]

Family Law Solutions

Prove Sobriety—Maintain Custody & Visitation Rights

In divorce, custody, or visitation cases, the parent's ability to care for children in the home becomes a concern when alcohol issues are called into question. Judges, family law attorneys, alternate dispute arbitrators, or social service agencies may regularly require parents to take sobriety tests. A parent may also volunteer to validate sobriety. Using the right alcohol detection tool results in sound custody and visitation decisions. The BI SoberTechTM handheld breathalyzer makes it easy to monitor sobriety, while driving compliance and accountability.

Benefits of Mobile Alcohol Detection

SoberTech is a lightweight, user-friendly alcohol detection device that discreetly prompts the user to test. This is convenient for individuals who may be at work or in a social setting. SoberTech uses biometric facial comparison technology to confirm client identity and the highest quality fuel-cell sensor to achieve an accuracy level within .005 Breath Alcohol Content (BrAC). Plus, near real-time delivery of results offers immediate peace of mind for families involved in custody or visitation arrangements.

Urinalysis results can take several days to process, and a negative result does not necessarily confirm abstinence, as alcohol can be quickly metabolized and eliminated from the body. With SoberTech, multiple tests can be conducted daily at scheduled or randomly times, based on the terms of the custody or visitation agreement. Test notifications are sent to the device, prompting the client to acknowledge receipt and submit a test. Test notifications can also be sent to the client's personal cell phone. Results are then conveniently sent to designated parties via text or email in real time.

KEY FEATURES

- · Discreet handheld device
- .02 BrAC or greater prompts a retest
- Scheduled or random testing
- Biometric facial comparison technology confirms client identity
- Professional grade fuel cell
- LCD color touch screen for easy reading and acknowledgment of test taking prompts
- Stores up to 500 tests
- Location point collected with time and date stamp

HOW IT WORKS

- Client enrolls, and testing schedule and alerts are set up
- Automated test reminders are sent to the device or a client's personal cell phone
- Client must test within 45 minutes of being prompted
- Near real-time results and alerts on missed or positive results are sent to designated recipients

For more information, call 800.876.1105 or visit us at bi.com

